INSTRUCTIONS FOR USE

Table of Contents

| | |
|--|-----------|
| <i>ENGLISH: INSTRUCTIONS FOR USE</i> | 2 |
| <i>SLOVENSKO: NAVODILA ZA UPORABO</i> | 8 |
| <i>ITALIANO: ISTRUZIONI L'USO</i> | 15 |
| <i>DEUTSCH: GEBRAUCHSANWEISUNG</i> | 22 |
| <i>HRVATSKO: UPUTSTVA ZA UPOTREBU</i> | 29 |
| <i>MAGYARORSZÁG: HASZNÁLATI UTASÍTÁS</i> | 36 |

ENGLISH: INSTRUCTIONS FOR USE

DEFINITION :

Mumu 3-ply surgical mask (facial mask for medical use), can be fit according to each face measures and shapes, flexible and can be used without disturbing the soft structure. Air permeable and lets breathing easily. Non-irritating. Provides protection against bacteria.

1. QUALITY SYSTEM OF MY TICARET VE MEDIKAL :

MY TICARET VE MEDIKAL has been manufacturing its products with the quality systems as given below;

- EN ISO 13485:2016 Quality Assurance System-Medical Devices

- EN ISO 9001:2015+AC Quality Assurance System

2. PRODUCT IDENTIFICATION AND RECOMMENDED USE :

The product is made by non woven fabric. The product's composition is polypropylene and does not includes latex. The product is breathable and has no special personal on environmental hazards. The product is made automatically in hygienic conditions. The product prevents the potential reactions between all kind of liquids and particles, microorganisms.

3. SUBSTANCE/MIXTURE OF RAW MATERIAL

| : Characteristic | Specification |
|------------------|--|
| Materials | |
| Mask Body | Outer Material Spunbond 25 gr (+/-5gr) |
| Filter Layer | Meltblown 25 gr (+/-5gr) |
| Inner Material | Spunbond 25 gr (+/-5gr) |
| Ear loop | Lycra-Elastane |
| Tie | Spunbond 40 gr (+/-5gr) |
| Nose-piece | Aluminium wire |

3.1. Material Safety Data Sheet (MSDS)

3.1.1. Composition:

Identification of the type of nonwoven product. Thermobonded Nonwoven by calendering process

Polypropylene CAS No: 25085-53-4 spunbond–meltblown modification

Web surface treatment - Concentration above 1°/0.....*No.

Binder.....*No

Additives.....*Yes.

May contain color pigments if the product is in color. White color may contain TiO₂. Other major components.....*No
Chemicals (in relevant concentration) that are in list of dangerous substances... *No

3.1.2. Hazards identification

Under normal conditions of use and handling, this product is not expected to create any health or safety hazards.

Accidental thermal decomposition or melting state can present hazards.

3.1.3. First-aid measures

Undernormal condition ;

Inhalation No specific measure to be taken
Skin contact No specific measure to be taken
Eyes contact..... No specific measure to be taken
Ingestion No specific measure to be taken

3.1.4. Fire fighting measures

1. Suitable extinguishing media;

Water spray, dry chemical or CO₂ extinguisher. No special procedures are expected to be necessary for this product. Normal fire fighting procedures should be followed to avoid inhalation of smoke and gases

2. Extinguishing media not to be used None

3. Special exposure hazards For flammable and toxic fumes as

well as skin contact with molten materials see § 10

4. Special protective clothing for fire-fighter..... None. It is recommended that fire-

fighters should wear full protective clothing including self contained breathing apparatus.

3.1.5. Accidental release measures

Personal Precautions: Avoid dust formation. Forms slippery surface.

3.1.6. Handling and storage

Keep in a dry and closed area with the original packing. The packages have to be handled so that they can not break and to be arranged as to prevent them from falling. The goods shall be handled with good industrial hygiene and safety practice.

3.1.7. Exposure controls / personal protection

No specific measures. Handle in accordance with good industrial hygiene and safety practice. Use of safety glasses and face mask is recommended if dust is formed during application.

3.1.8. Physical and chemical properties

Aspect.....Solid, in rolls or sheets

Appearance (the colour of the product as supplied).....Normally white if not a specific color is mentioned.

Odour.....Practically odorless

Ph.....Not applicable

Boiling point/boiling range Not applicable

| | |
|-----------------------------------|---|
| Melting point/melting range | (Polypropylene 165°C (330°F)) |
| Decomposition temperature | > 260°C (500°F) |
| Flash point | Not applicable |
| Flammability | Not easily flammable |
| Explosive properties | Not applicable |
| Oxidizing properties | Not applicable |
| Vapourpressure..... | Not applicable |
| Static electricity | The product can develop and/or accumulate static electricity, (i.e. by rubbing or friction) |
| Solubility | Mater insoluble - fat insoluble |
| Partition coefficient..... | Not applicable |

3.1.9. Stability and reactivity

The product is stable at room temperatures and does not decompose or self react when handled and stored under prescribed conditions. Toxic fumes can be generated under thermal decomposition.

3.1.10. Toxicological information

No toxic reaction known under normal conditions. Particularly, no case of coetaneous sensitisation or of mutagenic / carcinogenic activity is known. Under decomposition conditions, toxic fumes and contaminated water.

3.1.11. Ecological information

For transportation, storage and normal use no toxicological effect known. The fabric will not degrade biologically in short term.

3.1.12. Disposal considerations

As non hazardous solid waste, depending on local registration, nonwovens can be disposed of through recycling, landfill.

3.1.13. Transport information

Not classified as dangerous for transport.

3.1.14. Regulatory information

Not classified as dangerous in compliance with Turkey and European regulation regarding classifying, packaging and labeling of hazardous substances and products.

PERFORMANCE REQUIREMENTS OF FINISHED PRODUCTS

| We declare our products comply with EN 14683:2019+AC: 2019 as TYPE IIR and that all test results are the same for all the products (infant/adult/earloop/tie). Because all of the products are made from the same raw materials. | Units | Method | Results | Test Results |
|--|--------|-----------------------|---------|------------------------------|
| Parameters | | | | |
| BFE (Filtration) | % | EN 14683:2019+AC:2019 | ≥98 | % 98,4 (Ekoteks-20031582) |
| Differential | Pa/cm2 | EN 14683:2019+AC | <40 | 24,7 (Ekoteks-20031582) |

| | | | | |
|--|-------------|---|-----|---|
| pressure (Breathability) | | | | |
| Microbial cleanliness / | cfu/g | ISO 11737-1:2018 | ≤30 | 16 (Ekoteks-20031582) |
| bioburden (this is also requested in Portugal now) | | | | |
| Splash resistance pressure | kPa | ISO 22609 | ≥16 | 21,3 (Ekoteks-20031582) |
| Biocompatibility | ISO 10993-1 | Suitable for Skin Irritation, Cytotoxicity, Sensitization | | Cytotoxicity: It is not cytotoxicity. (Oxigen-2020-C-1099) Irritation: The sample has passed the irritation test. (Kocaeli University-2020-32) Sensitization: The sample has passed the sensitization test. (Kocaeli University-2020-33) |

5. INSTRUCTIONS FOR

USE 5.1. Intended Use:

Surgical mask (facial mask for medical use) is intended to be worn by medical personnel during surgical or other medical procedures to protect both the patient and the operating personnel and any other person that want and need to be protected, from transfer of microorganisms, body fluid, particulate material transfer and any other microbes.

Reduces exposure to blood and body fluids. Minimizes contamination to exhaled microorganisms. This product is intended for use in infection control practices.

5.2. Technical Specifications:

Non-irritating, Fluid Resistant, Three Ply construction. 3 pleats of folds to allow the user to expand the mask so it covers the area from the nose to the chin. Mask is secured with an ear loop to be placed behind the ears.

The surgical mask's (facial mask for medical use) three-ply layers work as follows:

- **The outer layer** repels water, blood, and other body fluids.
- **The middle layer** filters certain pathogens.
- **The inner layer** absorbs moisture and sweat from exhaled air.

5.3. Donning The Mask:

- Before putting on the mask, wash your hands for at least 20 seconds with soap and water, or rub your hands together thoroughly with alcohol-based hand sanitizer.
- Check for defects in the face mask, such as tears or broken loops.
- Position the colored side of the mask outward.
- If present, make sure the metallic strip is at the top of the mask and positioned against the bridge of your nose.
- Ear loops: Hold the mask by both ear loops and place one loop over each ear.

- Mold the bendable metallic upper strip to the shape of your nose by pinching and pressing down on it with your fingers.
- Pull the bottom of the mask over your mouth and chin.
- Be sure the mask fits snugly.
- Don't touch the mask once in position.
- If the mask gets soiled or damp, replace it with a new one.

!! Do not;

- touch the mask once it's secured on your face, as it might have pathogens on it
- dangle the mask from one ear
- hang the mask around your neck
- reuse single-use masks.

If you have to touch the face mask while you're wearing it, wash your hands first. Be sure to also wash your hands afterward, or use hand sanitizer.

5.4. Doffing The Mask:

- Before you take off the mask, wash your hands well or use hand sanitizer.
- Avoid touching the mask itself, as it could be contaminated. Hold it by the loops, ties, or bands only.
- Carefully remove the mask from your face once you: unhook both ear loops.
- Holding the mask loops discard the mask by placing it in a covered trash bin.
- After removing the mask, wash your hands thoroughly or use hand sanitizer.

FIRE FIGHTING MEASURES :

Suitable Fire Extinguishers and Methods:

Water spray, foam, carbon dioxide or dry chemicals. A sudden intervention should be made to the fire exit without any possible danger. If the material is melted, do not apply direct water flow. Use fine water spray or foam.

7. DISPOSAL CONSIDERATIONS :


Dispose of according to the Regulation on Control of Hazardous Wastes.

8. USE AREAS :

Hospitals, medical companies, doctor offices, laboratories, food manufacturers, cleaning companies and work places where the hygienic areas are necessary.

9. SHELF LIFE : 5 Years

10. SYMBOLS :

| | | | |
|---|------------------|---|---|
|  | Production Date |  | Shows the conformity to the European standards for the self-declared Class I products |
|  | Manufacturer |  | Keep Dry |
|  | Expiration Date |  | Keep out of sunlight |
|  | Lot Number |  | Non-Sterile |
|  | Reference Number |  | Do not use again |

11. STORAGE :

Capable of being stored continuously in ambient temperature of 10 to 30 deg C and relative humidity of 15 to 55%. Capable of operating continuously in ambient temperature of 10 to 30 deg C and relative humidity of 15 to 55%.

It is recommended that surgical face masks (facial mask for medical use)should be stored in their original containers and should be stored away from direct sunlight, heat sources and liquids, including chemicals. The area should be clean and should protect the masks from contamination. Never store it in a purse or pocket.

12. SAFETY INFORMATION :

- 📄 Pay attention to the warnings.
- 📄 It is not sterile.
- 📄 It is for single use only.
- 📄 Do not use if the package is damaged.
- 📄 Do not use the product after expiration date.

Declaration of conformity available on: <https://nundinae-group.com>

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